

Recombinant Anti-Human CD28 mAb

Catalog Number: C02302

Characteristics:

Source	Chinese Hamster Ovary cell line
Description	Isotype: Mouse IgG1/Kappa
Predicted molecular mass	150 kDa

Specification:


Appearance	Colorless clear liquid
Purity	≥95%, by SDS-PAGE (under reducing (R) & Non-reducing conditions, visualized by Coomassie staining)
Endotoxin	≤1 EU/mg by the LAL method
Activity	Measure its binding ability with Jurkat cells using flow cytometry. The ED ₅₀ of this effect is less than 30 ng/mL.
Formulation	0.22 μm-filtered solution containing PBS, pH 7.4
Size	1 mg/vial

Handling and Storage:


Shipping	The product is shipped with blue ice
Storage & Stability	24 months at -20°C to -80°C under sterile conditions
	6 months at 2°C to 8°C under sterile conditions
	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Data:

SDS-PAGE



Bioactivity



Recombinant Anti-Human CD28 mAb (Cat. No. C02302)
SDS-PAGE under reducing (R) & Non-reducing
conditions. The gel was stained with Coomassie Brilliant Blue G250.

Recombinant Anti-Human CD28 mAb (Cat. No. C02302)
binding with Jurkat cells. The ED₅₀ for this effect is less
than 30 ng/mL.

SEC-HPLC

