

Recombinant Human LIF

Catalog Number: Y03001

Characteristics:

Source Chinese Hamster Ovary cell line

Description Human LIF (Ser 23-Phe 202)
Accession # P15018

Predicted molecular mass 19.7 kDa

Specification:

Appearance White powder, Colorless clear liquid after reconstitution

Purity ≥95%, by SDS-PAGE (under reducing (R) & Non-reducing conditions, visualized by Coomassie staining)

Formulation Lyophilized from a 0.22 µm-filtered solution containing PBS, 5% mannitol and 0.01% Tween 80, pH 7.4

Size 50 µg/vial, 1 mg/vial

Handling and Storage:

Reconstitution It is recommended to redissolve in sterile deionized water.


Shipping The product is shipped with blue ice.

Storage & Stability 36 months at -20°C to -80°C in lyophilized state
6 months at -20°C to -80°C under sterile conditions after reconstitution
7-10 days at 2°C to 8°C under sterile conditions after reconstitution

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Data:

Bioactivity



Measured in a cell proliferation assay using TF-1 cells