

Recombinant Human LR3 IGF-1

Catalog Number: Y03101N

Characteristics:

Source	Chinese Hamster Ovary cell line
Description	Human LR3 IGF-1 (MFPAMPLSSLFVN+Gly49-Ala118(E51R)) Accession # P05019
Predicted molecular mass	9 kDa

Specification:


Appearance	White powder, Colorless clear liquid after reconstitution
Purity	≥95%, by SDS-PAGE (under reducing (R) & Non-reducing conditions, visualized by Coomassie staining)
Endotoxin	≤10 EU/mg by the LAL method
Activity	Measured in a serum-free cell proliferation assay using MCF-7 human breast cancer cells. The ED ₅₀ for this effect ≤1.5 ng/mL.
Formulation	Lyophilized from a 0.22 μm filtered solution in Acetic Acid.
Size	50 μg/vial, 1 mg/vial

Handling and Storage:

Reconstitution	It is recommended to redissolve in sterile deionized water.
Shipping	The product is shipped with blue ice.
Storage & Stability	36 months at -20°C to -80°C in lyophilized state 6 months at -20°C to -80°C under sterile conditions after reconstitution 7-10 days at 2°C to 8°C under sterile conditions after reconstitution Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Data:

Bioactivity



Measured in a serum-free cell proliferation assay using MCF-7 human breast cancer cells. The ED₅₀ for this effect ≤ 1.5 ng/mL. The activity of Y03101N is higher than other competing products.