

Recombinant Human BMP-2

Catalog Number: Y00421H

Characteristics:

Source	<i>E.coli</i>
Description	Human BMP2 (Gln283-Arg396) Accession # P12643
Predicted molecular mass	12.9kDa (monomer)

Specification:

Appearance	White powder, Colorless clear liquid after reconstitution
Purity	≥95%, by SDS-PAGE (under reducing (R) & Non-reducing conditions, visualized by Coomassie staining)
Endotoxin	≤100 EU/mg by the LAL method
Activity	Measured alkaline phosphatase production by ATDC-5 chondrogenic cells. The ED ₅₀ is ≤0.5μg/ml
Formulation	Lyophilized from a 0.22 μm filtered solution containing 0.1M Acetic acid, 5% Trehalose, pH 3.0
Size	10 μg/vial, 50 μg/vial, 1 mg/vial


Handling and Storage:

Reconstitution	It is recommended to redissolve in 4 mM HCl.
Shipping	The product is shipped with blue ice
Storage & Stability	36 months at -20°C to -80°C in lyophilized state 6 months at -20°C to -80°C under sterile conditions after reconstitution 7-10 days at 2°C to 8°C under sterile conditions after reconstitution Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Data:

V2.0

Bioactivity



Measured alkaline phosphatase production by ATDC-5 chondrogenic cells.