

Recombinant Human IL-17A

Catalog Number: Y06601H

Characteristics:

Source	Chinese Hamster Ovary cell line
Description	Human IL-17A(Ile20-Ala155) Accession #Q16552
Predicted molecular mass	15.5kDa

Specification:

Appearance	White powder,colorless clear liquid after reconstitution
Purity	≥95%, by SDS-PAGE (under reducing (R)& Non-reducing conditions, visualized by Coomassie staining)
Endotoxin	≤10 EU/mg by the LAL method
Activity	Measured by its ability to induce IL-6 secretion by NIH-3T3 mouse embryonic fibroblast cells. The ED ₅₀ for this effect is < 10 ng/mL.
Formulation	Lyophilized from a 0.22 µm-filtered solution containing PBS,5% Mannitol,0.01% Tween 80, pH7.4
Size	10 µg/vial, 50 µg/vial, 1 mg/vial

Handling and Storage:

Reconstitution	It is recommended to redissolve in sterile deionized water.
Shipping	The product is shipped with blue ice
Storage & Stability	36 months at -20°C to -80°C in lyophilized state 6 months at -20°C to -80°C under sterile conditions after reconstitution 7-10 days at 2°C to 8°C under sterile conditions after reconstitution Use a manual defrost freezer and avoid repeated freeze-thaw cycles.