

Recombinant Human IGF-1

Catalog Number: Y05301H

Characteristics:

Source	Chinese Hamster Ovary
Description	Human IGF-1 (Gly49-Ala118) Accession # P05019
Predicted molecular mass	7.7kDa

Specification:

Appearance	White powder, Colorless clear liquid after reconstitution
Purity	≥95%, by SDS-PAGE (under conditions, visualized by Coomassie staining)
Endotoxin	≤10 EU/mg by the LAL method
Activity	Measured in a serum-free cell proliferation assay using MCF-7 human breast cancer cells. The ED ₅₀ for this effect ≤60 ng/mL.
Formulation	Lyophilized from a 0.22 μm filtered solution containing PBS, 5% Mannitol and 0.01% Tween 80, pH7.4
Size	50 μg/vial, 1 mg/vial

Handling and Storage:

Reconstitution	It is recommended to redissolve in sterile deionized water.
Shipping	The product is shipped with blue ice
Storage & Stability	36 months at -20°C to -80°C in lyophilized state 6 months at -20°C to -80°C under sterile conditions after reconstitution 7-10 days at 2°C to 8°C under sterile conditions after reconstitution Use a manual defrost freezer and avoid repeated freeze-thaw cycles.